Product Certificate

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Qualifications & Patent

We have patent for LED mask , access to qualified EU market, advanced production systems, improved the quality of management

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Product quality inspection report

Our LED mask pass the test by China Quality Supervision and Inspection Center, Product safety non-toxic meet, safe household products, and skin and hair care appliances - Particular requirements, exceeding the national standard product requirements .

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Method one: Connect to the computer.
Take out the accessory cable of power supply, insert MICRO connector into DC-IN interface of power supply, and then connect the other end USB with the computer USB.
Method two: USB adaptor
Take out the accessory cable of power supply, insert MICRO connector into DC-IN interface of power supply, and then connect the other end USB with USB of DC 5V power supply.

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Matters needing attention

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Notices

- 1.Pregnant and minors are not suggested to use.
- 2.Part wavelength and frequency may cause epileptic disease, please do not use.
- 3.If suffering from thyropathy or taking thyropathy related medicines, please do not use.
- 4.If allergic to light or taking light-allergic medicines, for example,tetracyline, please do not use.
- 5.Try to avoid looking in the eye of the lighting area when use.

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Mask Maintenance

- 1.Use clear dry soft cotton cloth to wipe the surface gently in daily maintenance.
- 2.Do not use any solvent to clean the mask.
- 3.Do not disassemble the mask.
- 4.Please do not crush, throw, knock the mask and the power supply, keep power supply and lighting part dry.

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