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Light and Healthy

Light is a prerequisite for all living healthy growth

It's not only plants who need light . all creatures can not leave light . As early as the 2nd century we had the records of sunlight to treat disease .(solar therapy: sunlight, enhanced physique, prevent rickets).Artificial light therapy began at the end of the 18th century, to the middle of 19th century, visible light, infrared, ultraviolet, such as quickly formed, then the clinical treatment of widely in various fields of application and development.

LED Beauty Care

Light is proved to have medical properties

On October 17, 1903, the world's broadsheet newspaper published the Nobel foundation award announcements, on the winners list named "Niels Ryberg Finsen", has attracted a lot of people from all over the world care about the award ceremony Finsen is the person to create good phototherapy law lupus, smallpox, some with are patient his method to cure the patient, scrambling to flock to Stockholm, eager to see the world-famous imperial doctor.

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Application of modern optical care

Optical care is inspired by the light in the application of NASA and the Navy, a survey found that light can enhance the energy level, so that patients feel full of energy. Light is not intrusive, and it's low light pulses, three times brighter than the sun. Light can enhance the potential energy of skin cells, blue light can play sterilization and anti-inflammatory action to improve the surface skin, which has good effect for the treatment of acne and rosacea. Red light can promote the growth of skin collagen to achieve intensive skin care.

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LED Mask

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Introduction

LED optical mask made of high purity, high power density of red light, blue light to illuminate skin, can change the cellular structure, killing bacteria, to provide a suitable environment for new cells, enhance new collagen's elastin, promote cell growth, Complex inflammatory acne skin aging, and alleviate the sun burns skin, and not do harming the skin, can whiten skin, improve skin elasticity.

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Product principle-Wide spectrum of strong pulse photon technology

Specific wavelength pulse light energy can be through the epidermis, selective effects on the skin pigment, open blood vessels, hair follicles, atrophy, decomposition and absorption. At the same time also can make the thickening collagen tissue, generally after 15days to achieve whitening and anti-wrinkle treatment. Photon tender skin belongs to Non-ablative treatment and compare tender skin to before (chemical peeling, skinned paste), it won't cause adverse effects and complications, safe and reliable.

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
Product Schematic Diagram

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Features:

1. Non-contact, non-invasive treatment, quick effect, short course of treatment, no side effects.
2. Adjustable sleep, anti-inflammatory sterilization, effectively reduce acne.
3. The intensity of light source output stability, long life of light source.
4. No fever, no ultraviolet, won't produce pigmentation.
5. Can effectively strengthen the skin collagen cells activity, helps smooth tiny wrinkles.
6. Fade age spots, reducing solar sex skin damage.
7. Reduce the rosacea erythema and capillaries expansion.
8. Easy and convenient to operation.

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Introduction

As a new light source, with the biological effect is known by people. It plays an active role in the medical-grade beauty center and in the burn and scald treated instruments.

The international Medical has verified that the specific wavelength light can enhance the vitality of skin cell tissue to improve the skin condition by overlaid irradiation. Adopting the photoradiation theory, LED Mask can make the user maintain the facial skin safely and easily with intimate design. It can promote the growth of skin collagen to achieve greater depth cosmetology of skin. It is the phototherapy you can enjoy any time as you want. In addition, comparing with

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Create a natural beauty era

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Juvenility

Enhance the generation of collagen, reduce microgroove and wrinkles, and lighten skin color.

Compactness

Tighten the pores, improve the overall skin brightness, make the skin smooth and supple.

Rejuvenation

Fade freckles and opcs-redness, reduce acne and inflammation, accelerate the improvement of skin and the recovery after sunburn.

Moisturizing

Promote the circulation of face blood, lock up the skin aqua.

Vitality

Provide the cell activation energy, accelerate skin metabolism rate, thus to make the skin tender like a new-born baby.

Reduce acne

Restrain the inflamed phenomenon, accelerate the improvement of defective skin and repair of burn skin, and firm skin.

the traditional facial mask, LED Mask is more convenient, safer, and can be reused and has the advantage of environmental protection and others. LED Mask delivers voluminous LED red light energy to skin during the use, to process the crease from greater depth. It can not only flat the skin, but also prevent the growth of intradermal indigenous bacteria. You only need wear LED Mask 25 minutes every day, then LED red light will improve your skin gradually.

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Product data 产品参数

ITEM	Red Mask 2	Red Mask 1
Rated Voltage	DC5V±0.5V	DC5V±0.5V
Rated Current	1000 MA	500 MA
Power	5W	2.5W
Luminous Flux	300LUX	150LUX
Wavelength	650nm	650nm
Irradiated Area	30*20cm	30*20cm
Portable power supply Voltage	5V	5V

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Work Principle

LED red light adopting 650nm light efficient emitting, wavelength, through in the quantum photons absorbed by mitochondria of human body, and leads to the enzyme catalysis, so that to active the activity of many enzymes such as catalase and superoxide dismutase etc. Enhance metabolism of cells, improve immunity. Can stimulate and promote fibroblast proliferation and collagen synthesis. Red second generation Mask has timing function, it can be shut off automatically after using 20 minutes.

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Introduction

Blue light wavelength has the function for inhibiting melatonin. Light can affect the body's hormone secretion, heart rate, alertness, sleep quality, body and gene expression, blue-ray elevated body temperature, speed up the heart rate and reduce the role of sleepiness.

Blue light can inhibit sleep related delta brainwave, alpha brain-wave increased significantly related to alertness. This means that you can use the short wavelength light blue as lethargy of control measures. LED Blue-ray can have antiseptic and anti-inflammatory effects to improve the surface of the skin, have very good effect to treat acne and rosacea. Blue-ray is specifically kill cause acne acne bacillus, just

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Create a natural beauty era

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Bacteriostat

Repairing of acne damage and inhibition of wound bacteria.

Smoothing

Reduce the odds of scar caused by acne wound

Anti anaphylaxis

Slow down the skin inflammation

Glossiness

Improve the glossiness of skin.

Health

Reduce somnolence, increase vitality.

Conditioning

Adjust the biological clock (sleep disorders, time difference, menstrual syndrome).

need exposure 15-25 minutes on the skin, only for acne, light is very gentle, sensitive skin can be done, without any side effects. For the red and blue light acne r side effects, there is no need to worry about it.

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Product data 产品参数

ITEM	Blue Mask
Rated Voltage	DC5V±0.5V
Rated Current	500 MA
Output Power	2.5W
Luminous Flux	300LUX
Wavelength	470nm
Irradiated Area	30*20cm
Portable power supply Voltage	5V

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Work Principle

LED blue light is applied high luminescent wavelength of 470 nm, a specific blue light, with speed up blood circulation, activate cells, stimulate the active role of cell excitability. Can promote the synthesis of collagen, strengthen the skin collagen and elastic fibers. Can convey LED blue light to the skin, very comfortable, no concern of claustrophobia. Wearing the LED mask 25 minutes every day, your skin can be improved by LED blue light gradually.

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Three color mask

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
The LED mask inventor, invented the LED mask again!

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Collection of red and blue mask all characteristics, the effect is more strong

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

Introduction

Three colored light mask choose narrow-spectrum light is emitted from cold light and won't produce high fever, not burn the skin, it can turn the light energy into energy in cells, is the one of the significant instrument most safe for treat acne and skin.

It shows that red and blue mixed irradiation can amount to 35.7%, far higher than that of pure blue light group (8.3%) ($P=0.008$).

Severe acne is moderate acne have more obvious reaction, and the side effect is small, use blue light on the inflammatory response have to produce effect, cooperate with red curative effect is superior to use alone, no obvious side effects. Conclusion: the combination of red and blue light treatment of acne can obtain ideal curative effect, suitable for light treatment, patients with severe acne.

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
Change of power supply

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Red Mask
External power supply
红光面膜-外置电源

RBP Mask
Built-in power supply
Increase the Touch switch
三色光面膜-内置电源-增加开关

Power supply 电源
Touch switch 开关

Change of the product

1. The change of the product efficacy, multi-usage, satisfy to meet customer demand.
2. The change of product power supply, unibody design, facilitate better use for consumers.
3. The change of the light film, more large area light panel, make skin more fully by light area.
4. Increase of the regulator, to facilitate consumer use of products.

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Controller introduction

For better use and convenience, we improved our products and added regulator device, which is based on the use time of products, set up 6 gears products function. Touch the right of the mask holder, switch between various function.

- A level: Turn off automatically when the red and blue bright after 20min.
- B level: Turn off automatically when the red bright after 20min.
- C level: Turn off automatically when the blue bright after 20min.
- D level: Red bright longer.
- E level: Blue bright longer.
- F level: Power off.

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Product Certificate

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Qualifications & Patent

We have patent for LED mask, access to qualified EU market, advanced production systems, improved the quality of management

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Product quality inspection report

Our LED mask pass the test by China Quality Supervision and Inspection Center. Product safety non-toxic meet, safe household products, and skin and hair care appliances - Particular requirements, exceeding the national standard product requirements.

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Matters needing attention

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Notices

1. Pregnant and minors are not suggested to use.
2. Part wavelength and frequency may cause epileptic disease, please do not use.
3. If suffering from thyroopathy or taking thyroopathy related medicines, please do not use.
4. If allergic to light or taking light-allergic medicines, for example, tetracycline, please do not use.
5. Try to avoid looking in the eye of the lighting area when use.

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Mask Maintenance


1. Use clear dry soft cotton cloth to wipe the surface gently in daily maintenance.
2. Do not use any solvent to clean the mask.
3. Do not disassemble the mask.
4. Please do not crush, throw, knock the mask and the power supply, keep power supply and lighting part dry.

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Methods of battery charging

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Method one: Connect to the computer.

Take out the accessory cable of power supply, insert MICRO connector into DC-IN interface of power supply, and then connect the other end USB with the computer USB.

Method two: USB adaptor

Take out the accessory cable of power supply, insert MICRO connector into DC-IN interface of power supply, and then connect the other end USB with USB of DC 5V power supply.

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